Rosidal® sys

Ulcus cruris Kompressionsystem

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en Instructions for use

Product description and performance characteristics
Compression system for the treatment of venous leg ulcers over a period of approx. 12 weeks.

For an effective compression bandage, apply 2 Rosidal® soft and one Rosidal® K (8 cm and 10 cm) each time. Replacement bandages are included in this set to allow washing after each dressing change. Use the consumables (tg®, Porofix® and Mollelast® haft) as needed (see bandaging instructions, items 1, 9, 16 and 17). If these are used up more quickly than expected (e.g. in cases of extensive oedema), the components can also be bought separately.

Composition of the product

- 1x tg® (tubular bandage), size 7, 20 m, 67% cotton, 33% viscose (Ref. 24005)
- 4x Rosidal® soft (foam bandage, white), 10 cm x 2 m x 0.2 cm, 100% polyurethane (Ref. 13016)
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO), 8 cm x 5 m, 100% cotton (Ref. 22201)
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO), 10 cm x 5 m, 100% cotton (Ref. 22202)

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- Cut the tubular bandage to approx. 2.5 times the length from the end of the toes to the heel to keep and pull half of the bandage over the lower leg to below the knee. Have the patient hold the half that is still hanging off to keep the foot at an angle of 90 degrees to the ankle.
- Start applying Rosidal® soft from the top of the foot going around the inside of the instep including the base of the toes.
- The turns of the bandage should always overlap by approx. 1/2 - 2/3.
- Spiral wrap the heel including the heel.
- Return down around the ankle to secure the heel (towards the sole).
- To finish securing the heel, wrap above the ankle (towards the Achilles tendon).
- Continue to spiral wrap to the base of the calf muscles.
- Wrap the calf using a figure of eight pattern with Rosidal® soft.
- Secure the end of the bandage below the kneecap by wrapping the Mollelast® haft once around.
- Start applying the compression bandage Rosidal® K (8 cm wide) at the top of the foot going around the inside of the instep including the base of the toes.
- The turns of the bandage always overlap by approx. 1/2 - 2/3. Always apply Rosidal® K with full tension (max. stretch).
- Spiral wrap around the heel including the heel (towards the sole).
- To finish securing the heel, wrap above the ankle (towards the Achilles tendon) and continue to spiral wrap to the base of the calf muscle.
- Wrap the calf using a figure of eight pattern with Rosidal® soft.
- Secure the end of the bandage below the kneecap with a strip of Porofix®. Pull the tg® tubular bandage ending at the toes over the ankle to the knee.
- Secure the tg® tubular bandage below the knee and on the forefoot by wrapping Mollelast® haft around the ankle (towards the heel).
- The finished bandage.

Intended purpose

For venous leg ulcer compression therapy (see illustrations).

Indications

Venous leg ulcers and their comorbidities, as well as chronic oedema (due to venous insufficiency)

Contraindications

Advanced peripheral arterial occlusive disease, decompensated heart failure, septic phlebitis, plegmiasis coerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to product components.

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk.
- If the symptoms (pain and bluish tinge) do not disappear, the bandage must be removed and the patient informed by their doctor.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- With an Ankle Brachial Pressure Index (ABPI) greater than 0.8 or with an ABPI of 0.5-0.8 – use only under medical supervision.
- Always apply Rosidal® K with full tension (max. stretch).
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
- The width of the bandage should correspond to the diameter of the extremity to be banded.
- The bandage should be applied with tension around the leg in such a way that the pressure from the bandage decreases from the foot to the return.

Contraindications

Advanced peripheral arterial occlusive disease, decompensated heart failure, septic phlebitis, plegmiasis coerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to product components.

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk.
- If the symptoms (pain and bluish tinge) do not disappear, the bandage must be removed and the patient informed by their doctor.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- With an Ankle Brachial Pressure Index (ABPI) greater than 0.8 or with an ABPI of 0.5-0.8 – use only under medical supervision.
- Always apply Rosidal® K with full tension (max. stretch).
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
- The width of the bandage should correspond to the diameter of the extremity to be banded.
- The bandage should be applied with tension around the leg in such a way that the pressure from the bandage decreases from the foot to the return.

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Product description and performance characteristics
Compression system for the treatment of venous leg ulcers over a period of approx. 12 weeks.

For an effective compression bandage, apply 2 Rosidal® soft and one Rosidal® K (8 cm and 10 cm) each time. Replacement bandages are included in this set to allow washing after each dressing change. Use the consumables (tg®, Porofix® and Mollelast® haft) as needed (see bandaging instructions, items 1, 9, 16 and 17). If these are used up more quickly than expected (e.g. in cases of extensive oedema), the components can also be bought separately.

Composition of the product

- 1x tg® (tubular bandage), size 7, 20 m, 67% cotton, 33% viscose (Ref. 24005)
- 4x Rosidal® soft (foam bandage, white), 10 cm x 2 m x 0.2 cm, 100% polyurethane (Ref. 13016)
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO), 8 cm x 5 m, 100% cotton (Ref. 22201)
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO), 10 cm x 5 m, 100% algodón (Ref. 22202)

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