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## UK INSTRUCTIONS FOR USE

### Why pelvic floor training at all?

The pelvic floor is one of the most important muscle groups in our body. The pelvic floor trainer acts like a hammock. The muscles and ligaments hold the organs such as the bladder, intestines, vagina and uterus together elastically and stably on this "hammock". If this tension is released, the "hammock" can slide down. A strong pelvic floor can have a positive effect on your well-being and can increase your sense of pleasure.

### Quality that provides security.

The LUBEXXX® pelvic floor trainers weigh between 58 g and 92 g. They are made of high-quality elastomer or quality silicone and are very comfortable to wear. Easy to insert thanks to its ergonomic shape. Easy and discreet to use. With return strap.

### How does the pelvic floor trainer work?

Insert the pelvic floor trainer like a standard tampon. There are two weights inside the pelvic floor trainer. These begin to vibrate as soon as you physically move. The vibrations are transmitted to the muscles. The muscles try to balance the vibration. The pelvic floor muscle is trained. So all you have to do after insertion is move. (The 92 g trainer only has one weight).

### Application – see illustrations.

FIG. 1: Apply LUBEXXX® lubricant to the trainer as an insertion aid (if necessary).

FIG. 2: Insert the trainer into the vagina like a standard tampon. Find the right position.

Tipp, it's easier standing up.

FIG. 3: Now move while holding the trainer.

### Frequency of training.

Daily 1x to 2x approx. 10 minutes. If your muscles are weak, it may be a little difficult to hold the trainer at the beginning. Just keep going, because training can make your muscles stronger. Regular training is important. With successful training, the weight can be increased (advanced set with 72 g and 92 g).

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**EN LUBEXXX® Original** is a quality lubricant from Germany. Recommended by gynecologists. The skin compatibility was tested as “Very Good” (testing 03/2024, certificate at [www.lubexxx.de](http://www.lubexxx.de)). The water-silicone emulsion combines the advantages of a pure water lubricant and those of a pure silicone lubricant.

- Long-lasting gliding, plentiful
- Pleasantly silky skin feeling
- Water soluble, non-sticky
- Skin compatibility tested with “Very Good”
- Hormone-free
- Oil and grease free
- Easy dosing
- Suitable for rubber and latex
- Made in Germany

#### Application:

Simply apply to the desired parts of the body and enjoy the pleasantly silky skin feeling. Use more as needed. Simply wash off with water. Keep out of eyes.

#### Annotation:

This lubricant is not a medical device and is not intended for any medical purpose (e.g. the treat-

ment or relief of illnesses). Furthermore, the lubricant is not an accessory to a medical device and it is not intended to be used together with a medical device.

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