Black

Istr. HYALOTEND 20mg/2ml
Code/Revision: 33000477/A
Version: 4
Date: 25/08/20
Dimensions: 375 x 260 mm
Technical spec.: 581

fidia
 Fidia Farmaceutici S.p.A. - www.fidiapharma.com
 Under no circumstances should this artwork be altered
 without prior permission from Fidia Farmaceutici

Hyaluronic acid sodium salt solution for peritendinous use pre-filled syringe

PRODUCT INFORMATION

Hyalotend is a viscous solution of hyaluronic acid sodium salt (500-730 kDa) obtained by bacterial fermentation, in buffered physiological saline, for peritendinous injection. The solution (20mg/2ml) is contained in pre-filled syringe and is provided sterile. Hyaluronic acid is one of the main components of synovial fluid, moreover hyaluronic acid is also produced in the normal tendon sheath, and is a main component of tendon extracellular matrix. Administration of exogenous hyaluronic acid sodium salt in the peritendinous space, thanks to its viscoelastic properties, permits to reduce the surface friction of tendons, increases the gliding ability, reduces pain and improves function in tendinopathy thus speeding up return to normal activities and sport.

COMPOSITION

Principal component: Hyaluronic acid sodium salt
 Other components: Sodium chloride, Disodium hydrogen phosphate dodecahydrate, Sodium dihydrogen phosphate dihydrate, Water for injection.

INDICATIONS

Hyalotend is indicated for the treatment of pain and stiffness in tendinopathy of the upper and lower limbs.

ADMINISTRATION

Product administration should be performed exclusively by qualified physicians. Perform the administration with an ultrasound probe to guide the injection, if necessary. All the rules regarding the asepsis and the injection technique should be followed. Disinfectants containing quaternary ammonium salts should not be concomitantly used because hyaluronic acid can precipitate in their presence. Use suitable sterile needle (21-25 G) according to medical opinion. Discard the syringe and needle after single use.

DOSAGE

Administer 3 peritendinous injections of the product at weekly intervals, as appropriate according to medical opinion. More than one tendon may be treated at the same time, if necessary. Each syringe is intended for single tendon injection. Repeated treatments may be administered as required.

CONTRAINDICATIONS

Hyalotend should not be administered to patients with known hypersensitivity to any component of the product. Peritendinous injections of Hyalotend are contraindicated in cases of infections or skin diseases in the area of the injection site. Hyalotend has not been tested in pregnant women or individuals less than 18 years of age and therefore is contraindicated in these patient populations.

WARNINGS and PRECAUTIONS

Do not use the product after the expiry date reported on the package. Do not use the product if the packaging is opened or damaged: the sterility of the solution is guaranteed as long as the primary package is closed and undamaged. The syringe is for single use, that means it is intended to be used for one injection only. The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely administered. If a syringe is retained for a subsequent injection there is a risk of contamination that could lead to illness, infection and/or serious injury to the patient. If this product is reprocessed and/or reused, Fidia Farmaceutici cannot guarantee performance, functionality, material structure, cleanliness or sterility of the product. Reprocessing and or reusing could cause serious harm to the patient health and safety. After use, dispose according to applicable national practice. As with any invasive joint procedure, it is recommended that care should be taken not to overburden the joint immediately after the peritendinous injection. Keep out of reach of children.

UNDESIRABLE EFFECTS

Local reactions such as pain, and warmth or redness/ swelling /bruising at the injection site might be expected. Usually such symptoms are transient and disappear spontaneously within a few days by resting the affected joint and applying ice locally. Local or systemic allergic reactions due to individual hypersensitivity have been sporadically recorded after hyaluronic acid injection.

STORAGE

The product must be stored at temperature not exceeding 25°C in the original package protected from light. Do not freeze.

HOW SUPPLIED

Box containing 1 pre-filled syringe sealed in blister.
 Box containing 3 pre-filled syringes, each sealed in blister.

The solution of 20 mg/2 ml hyaluronic acid sodium salt contained in the pre-filled syringe is sterilised using steam.

Manufacturer FIDIA Farmaceutici S.p.A.
 Via Ponte della Fabbrica 3/A
 35031 Abano Terme (PD), Italy

Date of the latest revision of the instructions for use: May 2018

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Výrobca FIDIA Farmaceutici S.p.A.