L-Mesitran® TULLE Non-adhering antibacterial dressing with L-Mesitran Soft

NEDERLANDS

L-Mesitran Tulle

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L-Mesitran Tulle

Product description

L-Mesitran Tulle is a non adherent polyethylene dressing im-pregnated with the patented L-Mesitran Soft gel that contains: 40% medical grade honey, medical grade hyaloïpergenic lanolin (Medilan®), propylene glycol, PEG 4000 and vitamins C & E. When the gel is in contact with the wound, fluid is drawn from the surrounding tissues. By doing so, a moist wound healing envi-ronment is created. This environment stimulates the wound healing process: by facilitating autolysis of necrotic and other devitalised material, aids wound bed preparation and promotes better epi-thelial cell migration.

- L-Mesitran Tulle has antibacterial properties.

- The gel will dilute gradually and can be removed easily as per local wound cleansing protocols.

- L-Mesitran Tulle is a primary wound dressing that should be covered with most commonly used secondary dressings.
- The properties of L-Mesitran Tulle prevent the secondary dress-ing adhering to the wound bed.
- Due to the use of natural ingredients, each batch of gel may vary in colour (white - yellow - brown).
- The gel may contain visible sugar crystals; however this will not affect the effectiveness and/or application of the gel.
- The white polyethylene film will allow exudates to pass to the secondary dressing.

Indications

Superficial wounds, acute wounds such as cuts, abrasions and donor sites. Burn or laser wounds, superficial and partial thickness burns (1st & 2nd degree). Chronic wounds such as pressure ulcers, venous, arterial and diabetic ulcers. Fungating wounds (to help deodorise and/or debride). Colonised acute wounds and (post-operative) surgical wounds.

Contraindications and constraints

No contraindications are known to date. Do not use on patients who are known to be sensitive to the dressing or any of its components.

Preventive measures and comments

- Do not use if the packaging is damaged or already opened. Do not use after the expiry date.
- Consult a doctor if further measures are deemed necessary for proper treatment of the wound.
- Inspect and clean the wound according to local protocols before use.
- Due to autolytic debridement the wound may appear to increase in size after the first dressing changes. This is normal and occurs as the devitalised tissue is removed, preparing the wound bed for healing.
- Seek medical supervision if the patient shows signs of infection (fever, pain) or if the condition of the wound deteriorates.
- Each L-Mesitran Tulle dressing is for single use only. Dispose the used product as per local protocol.

Instructions for use

- Inspect and if necessary, cleanse the wound as per local protocols.
- Remove L-Mesitran Tulle from the pouch and cut it to the de-sired size if necessary.
- Peel the dressing off both paper backing sheets and apply L-Mesitran Tulle to the wound.
- Cover and secure with the appropriate secondary dressing de-pending on the amount of wound fluids.
- Inspect the wound and the surrounding tissue regularly, as per local protocols.
- In the event of intolerance to the dressing or one of its compo-nents, remove the dressing using appropriate techniques and clean the area carefully.

Frequency of change

L-Mesitran Tulle can remain in place for several days, dependent upon the amount of exudate and the needs of the wound. Re-place L-Mesitran Tulle when the gel has gone. Where appropriate, L-Mesitran Tulle may remain in place while the secondary dress-ing is changed. In the event L-Mesitran Tulle has become too dry, complete and gentle removal of L-Mesitran Tulle may be aided by irrigation.

Storage precautions

Store at room temperature between 5°C and 25°C (41°-77°F). Avoid direct sunlight and moisture. L-Mesitran is a natural product and the gel may become liquid under higher temperatures. If this is expected to be a problem, then the dressing (without the outer carton) may be placed in a fridge prior to use.

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