Geratherm[®] VIVO

Arm Blood Pressure Monitor

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DE

Instructions for Use

EN

Mode d'emploi

FR

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IT

Manual de instrucciones

ES

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PT

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NL

DE SYMBOLERKLÄRUNG **EN SYMBOL DEFINITIONS** **FR DÉFINITIONS DES SYMBOLES** **IT DEFINIZIONI DEI SIMBOLI** **ES DEFINICIONES DE SÍMBOLOS**
PT DEFINIÇÕES DE SÍMBOLOS **NL SYMBOOLDEFINITIES**

	Warnhinweise, siehe beiliegendes Dokument, Warning information, refer to the attached document, Informations sur les avertissements, voir le document ci-joint, Per informazioni sulle avvertenze, consultare il documento allegato, Informação sobre advertências, consulte el documento adjunto, Para informações de aviso, consulte o documento anexo, Informatie over waarschuwingen, zie het bijgevoegde document
	Geräteklassifizierung: Anwendungsteil Typ BF, Device classification: type BF applied part, Classification du dispositif : type BF partie appliquée, Classificazione del dispositivo: tipo BF parte applicata, Clasificación del dispositivo: tipo BF parte aplicada, Classificação do dispositivo: peça tipo BF aplicada, Apparaat classificatie: toegepast deel type BF
	Gerät darf nicht im Hausmüll entsorgt werden. The device must not be disposed of with household waste. Ne pas éliminer l'appareil avec les ordures ménagères. Non smaltire il dispositivo nei rifiuti domestici. El aparato no debe desecharse con la basura doméstica. O aparelho não deve ser eliminado juntamente com o lixo doméstico. Apparaat mag niet bij het restafval.
	Gebrauchsanweisung beachten, Follow the instructions for use, Suivre la notice d'utilisation, Seguire le istruzioni per l'uso, Observe las instrucciones del manual de uso, Consulte as instruções de utilização, Raadpleeg de gebruiksaanwijzing
	Trocken aufbewahren, Keep dry, A protéger contre les projections d'eau, Proteggere dall'umidità, Proteger contra la humedad, Manter seco, Droog bewaren
	Vor Sonnenlicht schützen, Keep away from the sunlight, Protéger contre la lumière du soleil, Proteggere dalla luce del sole, Proteger de la luz solar, Manter afastado da luz solar, Uit de buurt van zonlicht houden
IP21	Das Gerät ist gegen Spritzwasser geschützt. Wasser, das aus jeder Richtung gegen das Gehäuse spritzt, darf keine schädlichen Wirkungen haben. The device is protected against splashing water. Water splashed against the enclosure from any direction shall have no harmful effects. L'appareil est protégé contre les éclaboussures d'eau. Le boîtier dans n'importe quelle direction ne doit pas avoir d'effets nocifs. L'apparecchio è protetto contro gli spruzzi d'acqua. L'acqua spruzzata contro l'involucro da qualsiasi direzione non avrà effetti dannosi. El aparato está protegido contra salpicaduras de agua. El agua que salpique contra el recinto desde cualquier dirección no tendrá efectos nocivos. O dispositivo está protegido contra salpicos de água. Os salpicos de água no compartimento a partir de qualquer direção não têm efeitos prejudiciais. Het apparaat is beschermd tegen opspattend water. Water dat vanuit welke richting dan ook tegen de behuizing spat, zal geen schadelijke gevolgen hebben.
	Hersteller, Manufacturer, Fabricant, Produttore, Fabricante, Fabricante, Fabrikant
	Herstellungsdatum, Date of manufacture, Date de fabrication, Data di produzione, Fecha de fabricación, Data de fabrico, Datum van vervaardiging

DE SYMBOLERKLÄRUNG EN SYMBOL DEFINITIONS FR DÉFINITIONS DES SYMBOLES IT DEFINIZIONI DEI SIMBOLI ES DEFINICIONES DE SÍMBOLOS PT DEFINIÇÕES DE SÍMBOLOS NL SYMBOOLDEFINITIES

LOT	Chargenbezeichnung, Batch code, Désignation du lot, Designazione del lotto, Denominación de lote, Código do lote, Batchcode
EU REP	Bevollmächtigter Vertreter in der Europäischen Gemeinschaft/ Europäischen Union, Authorized representative in the European Community/ European Union, Représentant autorisé dans la Communauté Européenne/l'Union Européenne, Rappresentante autorizzato nella Comunità Europea/Unione Europea, Representante autorizado en la Comunidad Europea/Unión Europea, Representante autorizado na Comunidade Europeia/União Europeia, Gemachtigde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie
CH REP	Bevollmächtigter Vertreter in der Schweiz, Authorized representative in Switzerland, Représentant autorisé en Suisse, Rappresentante autorizzato in Svizzera, Representante autorizado en Suiza, Representante autorizado na Suíça, Geautoriseerde vertegenwoordiger in Zwitserland
UK REP	Bevollmächtigter Vertreter im Vereinigten Königreich, Authorized representative in United Kingdom, Représentant autorisé au Royaume-Uni, Rappresentante autorizzato nel Regno Unito, Representante autorizado en el Reino Unido, Representante autorizado no Reino Unido, Gemachtigde vertegenwoordiger in het Verenigd Koninkrijk
MD	Medizinprodukt, Medical Device, Dispositif médical, Dispositivo medico, Dispositivo médico, Dispositivo médico, Medisch apparaat
#	Modellnummer, Model number, Numéro de modèle, Numero modello, Número de modelo, Número do modelo, Modelnummer
UDI	Eindeutige Produktidentifizierung, Unique device identification, Identifiant unique de l'appareil, Identificatore univoco del dispositivo, Identificador único de dispositivo, Identificador único do dispositivo, Unieke apparaatidentificatie
	Lagertemperatur zwischen -20 °C und 55 °C, Storage temperature between -20 °C and 55 °C, Température de stockage entre -20 °C et 55 °C, Temperatura di stoccaggio compresa tra -20 °C e 55 °C, Temperatura de almacenamiento entre -20 °C y 55 °C, Temperatura de armazenamento entre -20 °C e 55 °C, Bewaar temperatuur tussen -20 °C en 55 °C
	Lagerung bei einer relativen Luftfeuchtigkeit zwischen 10 % und 93 % R.F. Storage at a relative humidity between 10 % and 93 % RH. Stockage à une humidité relative de l'air entre 10 % et 93 % RH. Conservazione ad un'umidità relativa compresa tra il 10 % e il 93 % RH. Almacenamiento a una humedad relativa entre 10 % y 93 % RH. Armazenamento a humidade relativa entre 10 % e 93 % de HR. Opslag bij een relatieve luchtvochtigheid tussen 10 % en 93 % RH.

DE SYMBOLERKLÄRUNG EN SYMBOL DEFINITIONS FR DÉFINITIONS DES SYMBOLES
 IT DEFINIZIONI DEI SIMBOLI ES DEFINICIONES DE SÍMBOLOS
 PT DEFINIÇÕES DE SÍMBOLOS NL SYMBOOLDEFINITIES

	Lagerung bei einem atmosphärischen Druckbereich zwischen 70 kPa und 106 kPa. Storage at an atmospheric pressure range between 70 kPa and 106 kPa. Stockage dans une plage de pression atmosphérique de 70 kPa à 106 kPa. Stoccaggio ad un intervallo di pressione atmosferica compreso tra 70 kPa e 106 kPa. Almacenamiento a un rango de presión atmosférica entre 70 kPa y 106 kPa. Armazenamento em intervalo de pressão atmosférica entre 70 kPa e 106 kPa. Opslag bij een atmosferisch drukk bereik tussen 70 kPa en 106 kPa.
	Importeur, Importer, Importateur, Importatore, Importador, Importador, Importeur
	Gleichstrom, Direct Current, Courant direct, Corrente continua, Corriente continua, Corrente contínua, Gelijkstroom

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CE 0123

TABLE OF CONTENTS

General Description	23
Warranty	23
Unpacking Inspection	23
Package Content	23
Intended Use	23
Contraindications	23
Warnings and Notes	24
Description of the device	26
WHO Blood Pressure Indicator	27
Power Connection	27
Function Setting	28
How to Apply the Cuff	29
How to Take Proper Measurements	30
Error messages and Handling	32
Cleaning and Disinfection	34
Questions and Answers on Blood Pressure	35
Specifications	36
Authorized Component	37
Appendix: EMC Information	37

GENERAL DESCRIPTION

Thank you for purchasing the Arm Blood Pressure Monitor Geratherm® vivo. Designed for ease of use and accuracy, the device detects your blood movement through your brachial artery and converts the movements into a digital reading based on oscillometric method to measure blood pressure. We recommend that you read through this manual carefully before using the device for the first time.

WARRANTY

The device is under warranty for 24 months from the date of purchase. The batteries, the packaging, cuff and any damage caused by improper use are not covered by the warranty. Excluding the following user-caused failures:

- 1) Failure resulting from unauthorized disassembly and modification.
- 2) Failure resulting from an unexpected dropping during application or transportation.
- 3) Failure resulting from not following the instructions in the operating manual.

UNPACKING INSPECTION

Before using, please open the package carefully and check if all the parts are available according to the following packing list and if there are any damaged parts during transportation, then install and operate strictly in accordance with the user manual.

PACKAGE CONTENT

- 1x Blood Pressure Monitor
- 1x Cuff 22 - 42 cm (8.6" - 16.5")
- 4x Alkaline batteries type AAA 1.5 V
- 1x Storage Bag
- 1x Instructions for Use
- 1x USB-Type-C Cable

EN

INTENDED USE

The Arm Blood Pressure Monitor is intended to measure the systolic blood pressure and diastolic blood pressure, as well as the pulse rate of the adult person via non-invasive oscillometric technique at medical facilities or at home.

CONTRAINDICATIONS

This device is not for infants or mentally-impaired individuals who cannot express their thoughts.

WARNINGS AND NOTES

Warnings:

- Do not use this monitor for infants, toddlers, children, or persons who cannot express themselves.
- Do not use this monitor for any purpose other than measuring blood pressure and pulse rate.
- Do not take or make adjustments to any medications by yourself based on readings from the device. Contact your physician for specific information about your blood pressure. The patient should not self-diagnose or self-medicate per measured results. Kindly adhere to the instructions of your physician or health provider.
- Do not apply the cuff on an arm with an unhealed wound.
- Do not apply the cuff on an arm receiving an intravenous drip or blood transfusion.
- Do not disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.
- The device and cuff are not water resistant. Prevent rain, sweat and water from soiling the device and cuff.
- Do not use this monitor in a moving vehicle or other unstable environment.
- Precision components are used in the construction of this device. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.
- Do not use this device adjacent to or stacked on other jamming equipment to prevent improper functioning. If such use is necessary, all accessories involved must be checked for normal operation.
- Do not use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- To obtain the greatest accuracy from your blood pressure device, it is recommended that the device be used within the specified temperature and the relative humidity, please see the Technical Specifications.
- The use of other accessories other than those specified or provided by the equipment manufacturer may cause electromagnetic radiation to increase or decrease electromagnetic immunity resulting in operational failure.
- Do not attempt to service or repair this device yourself. Should a malfunction occur, refer to local distributor or the manufacturer.


WARNINGS AND NOTES

Notes:

- This device shall not be serviced or maintained when using.
- Maintenance should be done by the manufacturer as suggested.
- NEVER diagnose or treat yourself based on the readings. ALWAYS check with your doctor about any possible causes and explanation about blood pressure values. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.
- Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pregnancy, pre-eclampsia, or renal disease.
- Consult with your physician before using this monitor if you have had a mastectomy.
- Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders, because the cuff inflation can cause bruising.
- The pulse display is not suitable for checking the frequency of heart pacemakers! In cases of irregular heartbeat, measurements made with this device should only be evaluated after consultation with your doctor.
- Avoid tightly folding the cuff or storing the air hose tightly twisted for long periods, as such treatment may shorten the life of the components.
- To avoid strangulation, please keep the air tube and power supply line away from the infants, toddlers, children, and pets.
- Do not use the batteries and the AC adapter to provide power at the same time. Remove the battery if the device is not likely to be used for some time.
- Used equipment, parts and batteries are not treated as ordinary household waste and must be disposed of according to the applicable local regulations.


DESCRIPTION OF THE DEVICE

Blood pressure monitor


This product is composed of the blood pressure monitor and cuff.

Display Screen


1. Date and Time
2. Systolic Blood Pressure Icon
3. Systolic Blood Pressure Value
4. Diastolic Blood Pressure Icon
5. Diastolic Blood Pressure Value
6. Heartbeat Icon
7. Irregular Heartbeat Icon
8. Pulse Rate Icon
9. Memory Icon
10. Pulse Rate Value
11. Unit of Blood Pressure
12. Correct Cuff Wrapping
13. Blood Pressure Indicator
14. "Keep Still" Indication
15. User Icon
16. Guest
17. Battery Icon
18. Average Value Icon


WHO BLOOD PRESSURE INDICATOR


White
for Standby

Green
Bar Graph for Normal
Blood Pressure

Yellow
Bar Graph for
Mild-high/ Mild-low
Blood Pressure

Red
Bar Graph for
High Blood Pressure

Systolic Blood Pressure (mmHg)	Hierarchical Relationship	Diastolic Blood Pressure (mmHg)	Color Display
≥ 160	or	≥ 100	red
140 - 159	or	90 - 99	yellow
90 - 139	or	60 - 89	green
< 90	and	< 60	yellow


Warning: Never diagnose or treat yourself based on the readings. Please always consult with your physician.

POWER CONNECTION

1. Install the Batteries

- Open the battery cover according to the method shown in the figure.
- Place 4 AAA batteries 1.5 V in the battery compartment and pay attention to the electrode indication of the batteries. Install the battery as indicated in the picture.

Remove the batteries if you do not intend to use the device for a long time (over 3 months).

EN

2. Type-C Connection for Power Supply

In addition to the batteries, the power can also be supplied by plugging into 5 V external power supply through a Type-C port. (Charging Adaptor is not included in the package).

External Power
Supply Port

FUNCTION SETTING


1. User Setting

a) Guest Mode

In the power-off mode, press the "⌘" button to enter Guest Mode "👤". Under the guest mode, the blood pressure readings will not be recorded in the device.

b) User 1/ User 2 Mode

Then press the "⌘" button again to switch and select User 1 "👤" or User 2 "👤". Under the User 1/2 Mode, the blood pressure readings will be recorded in the device with date and time. Be certain that the same user is selected for specific person for continuous reading memory and easy checking.


2. Date and Time Setting

In the power-off mode, press the "⌘" button for about 3 seconds to enter the date setting interface, and the "year" will flash.

Press the "🔍" button to adjust the year and press the "⌘" button to confirm the selection.


When the year is set, it will automatically enter the month setting. Press the "🔍" button to adjust the month and press the "⌘" button to confirm the selection. Follow the same steps to date/hour/minute.

FUNCTION SETTING


3. Unit Display Setting

There are two units of blood pressure display, mmHg and kPa. The default unit is mmHg. In the power-off mode, press the "START/STOP" button for about 5 seconds to enter the unit selection. Press the "←" button to switch between mmHg and kPa, and then press the "OK" button to confirm the selection.


HOW TO APPLY THE CUFF

1. Connect the cuff to the monitor by inserting the air plug into the air jack.

2. Place your hand through the cuff loop. Pull the cuff until it reaches your upper arm.


Note:

The bottom edge of the arm cuff should be 2-3 cm (0.8" - 1.2") above the elbow. The air tube should be on the inside of your arm and aligned with your middle finger.

Make sure that the air tube is positioned on the inside of your arm and wrap the cuff securely, so it cannot move around your arm.

Note:

Repeated measurement will result in blood congestion in the arm, which will affect the measurement result.

HOW TO APPLY THE CUFF

How to avoid blood congestion and ensure the repeated measurement is accurate?

You can raise the left hand and hold the fist several times or take off the cuff and rest for at least 5 minutes before taking the measurement.

3. Sitting correctly

To take a measurement, you need to be relaxed and comfortably seated in a room with a comfortable temperature.

- Sit in a comfortable chair with your back and arm supported.
- Keep your feet flat and your legs uncrossed.
- The arm cuff should be placed on your arm at the same level as your heart, with the arm resting comfortably on a table.

HOW TO TAKE PROPER MEASUREMENTS


1. Preparation before Measurement

- Take off the clothes on the arm.
- Always measure in the same arm (generally the left arm).
- Remain still during the measurement.
- Relax as much as possible and do not talk during measurement.
- Measure your blood pressure at about the same time every day.
- Do not measure the blood pressure immediately after physical exercise or a bath. Rest for 20 - 30 minutes before taking the measurement.
- Measurements under the conditions listed below may affect results: eating, alcohol consumption; drinks containing caffeine, sport, talking, being nervous, being in unsteady mood, bending forward, moving, room temperature changes, being inside a moving vehicle, repeated and continuous measuring.

HOW TO TAKE PROPER MEASUREMENTS

2. Taking a Measurement

- 1) Fasten the arm cuff following the instruction of "How to Apply the Cuff". Start the measurement after wearing the cuff correctly.
- 2) Press the "START/STOP" button. The monitor will start inflating for measurement and display "00". Check the measured values after the measurement finished.


Note:

If you feel uncomfortable during the measurement, press the "START/STOP" button immediately to stop the measurement. Please consult your doctor if unexpected readings are obtained.

3. Memory Function


- 1) Each measured value is stored automatically under the appropriate user group. This device can store up to 199 sets of measurements for User 1 and User 2. Once the memory log is full, old values will be refreshed with new ones.
- 2) In the power-off mode, press the "MEMO" button once and the device will display the average value of the blood pressure measurements of the last 2 or 3 times. Press the "MEMO" button again, and the latest measured value will be displayed. Press the "MEMO" button again and the rest measured values will be displayed one by one.

HOW TO TAKE PROPER MEASUREMENTS

3) Delete Memory

- 1 In the power-off mode, press the "⌂" button to select the user group whose measured values need to be deleted.
- 2 Press the "START STOP" button to power off the device and press the "⌂" button once to activate the screen.
- 3 Keep pressing the "⌂" button for about 3 seconds to delete the memories of the selected user and the "R3" icon will appear on the screen.

4. "Cuff Wrapping" Detection

The "⌂" icon is displayed stationary on the screen when the cuff is wrapped correctly. When the cuff is too loose, the "⌂" icon will always flash to remind you. If the "⌂" icon is flashing, please press the "START STOP" button to stop the measurement.

5. "Keep Still" Indication

The "⌂" icon flashes when you move your body or shake your arm during the measurement, which may cause incorrect measurement results. Please adjust your posture and measure again.

ERROR MESSAGES AND HANDLING

If the measurement is faulty, any of the following error codes may appear. Kindly use the recommended method for measurement.

Errors	Cause	Solution
Er U	The pressure cannot reach 30 mmHg (4 kPa) in 12 seconds.	Please make sure that the cuff is inserted and worn correctly. Please check that the cuff and tube is not broken.
Er H	The inflation reaches 295mmHg, and it deflates automatically after 20ms.	Please make sure that the cuff is worn correctly, and the measuring posture is proper.
Er 1	The pulse rate is not detected correctly.	Please make sure that the cuff is worn correctly, and the measuring posture is proper. Please check that the measuring arm is free of clothing.
Er 2	Too much disturbance (moving, talking or magnetic disturbance during a measurement).	Please do not move or talk during measuring.

ERROR MESSAGES AND HANDLING

Er 3	The measurement result is abnormal.	Please do not exercise strenuously before measurement. Please do not move or talk during measuring.
Er 23	SYS value is lower than 57 mmHg.	Please make sure that the cuff is worn in the right position and correctly. Please do not move or talk during measuring.
Er 24	SYS value is higher than 255 mmHg.	
Er 25	DIA value is lower than 25 mmHg.	
Er 26	DIA value is higher than 195 mmHg.	

Troubleshooting

Problem	Possible Faulty	Solution
Failure to power on	Whether the power is insufficient.	Replace the batteries or insert the cable for power supply.
	Whether the positive and negative poles of the battery are installed reversely.	Install the batteries correctly.
No pressurizing	Whether the air tube plug is inserted tightly.	Insert the air tube plug firmly into the jack.
	Whether the air tube is broken or leaked.	Please contact the dealer to replace with a new cuff.
The reading is extremely low (or high)	Whether the arm is moved when pressurization.	Keep your arm and body still.
	Whether you talk during the measurement.	Keep quiet while measuring the blood pressure.
Air leakage of the cuff	Whether the cuff is twined too loose.	Please tighten the cuff.
	The airbag of the cuff is ripped.	Please contact the dealer to replace with a new cuff.
 If the blood pressure still cannot be measured after trying the above-stated solutions, please contact the dealer. Do NOT attempt to disassemble the device by yourself.		

CLEANING AND DISINFECTION

1. Cleaning

The device can be cleaned with a soft, clean cloth dampened with a small amount of neutral detergent or water.

- ⚠ Do not use corrosive cleaning agents and take care not to immerse any part of the monitor in any fluid.
- ⚠ Do not allow water or other liquids to flow into the device.

2. Disinfection

Recommended Disinfecting Agent

75 % medical alcohol

Carefully wipe the device with a soft, clean cloth dampened with a small amount of the above disinfectant, and dry immediately with a soft, clean, dry cloth.

- ⚠ Do not disinfect through methods like high-temperature steam or ultraviolet radiation. These might damage the device and reduce its service life.

It is suggested to disinfect the monitor before and after use each time. Each time of disinfection shall be completed within 1min. The number of repeated disinfections each time shall not exceed 2 times.

3. Disposal

Dispose of the monitor, other components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.

Notes

- Do not bend or crease the air tube excessively.
- Do not store the monitor or its components:
 - if the monitor or its parts is wet.
 - in locations with extreme temperatures, humidity, direct sunlight, dust, or corrosive gases.
 - in areas with a high risk of vibrations or shocks.

QUESTIONS AND ANSWERS ON BLOOD PRESSURE

Why is the blood pressure value obtained at home lower than that obtained at the hospital?

- The blood pressure difference between home and hospital measurements is about 20 mmHg - 30 mmHg (2.7 kPa - 4.0 kPa). This is because individuals tend to be more relaxed at home than at the hospital.
- In addition, when the device is placed at a position over the heart, the blood pressure value tends to be much lower than it actually is. Ensure the device is positioned right at the heart level.

Why is the blood pressure value obtained at home higher than that obtained at the hospital?

- The anti-hypertensive drug might have lost its efficacy. Kindly adhere to your doctor's instructions.
- The cuff might not be in the correct position. If the cuff is not placed right, no arterial pressure value will be obtained, and the blood pressure value might be much higher than it is. Therefore, properly position the cuff.
- The cuff is not tight enough. If the cuff is loose, the compression force might fail to transmit to the artery, causing the blood pressure value to be much higher than it is. Therefore, re-adjust and tighten the cuff further.
- The patient is not sitting correctly during the measurement. Slouching, tilting, bending, and sitting cross-legged are not encouraged while taking blood pressure measurements due to increased abdominal pressure or the arm position being below the heart. Kindly take readings in the correct posture.

When can I obtain better measurements?

- Measurements are best taken in the mornings right after you urinate or when your mind and body are stable. We recommend taking readings at the same time of the day, every time.

SPECIFICATIONS

Model	ARM-30T	
Display	LED screen	
Measuring Method	Oscillometric measurement	
Measuring Part	Upper arm	
Pneumatic Pressure Measuring Range	0 - 295 mmHg (0 - 39.3 kPa)	
Maximum Pressure Protection	295 mmHg (39.3 kPa)	
Measurement Range	Blood pressure value	SYS: 57 - 255 mmHg (7.6 - 33.4 kPa) DIA: 25 - 195 mmHg (3.33 - 26 kPa)
	Pulse rate	40 - 199 beats per minute
Accuracy	Blood pressure value	±3 mmHg (±0.4 kPa)
	Pulse rate	±5 %
Low Battery	When the power is lower than 4 V ±0.1 V, the device will be turned off automatically.	
Power Source	4 x AAA batteries or 5 V  Type-C cable	
Memory	2 users x 199 memories + guest mode	
Dimension	130 mm (L) x 100 mm (W) x 49 mm (H) (5.12" x 3.94" x 1.93")	
Screen Size	67 mm (L) x 68 mm (W) (3.7")	
Cuff Size	22 - 42 cm (8.6" - 16.5")	
Weight	About 221 g (without batteries)	
Auto Power-off	1 minute without operation	
Anti Electronic Shock Degree of Protection	Type BF	
Operation Mode	Continuous operation	
Protection Against Harmful Ingress of Water or Particular Matter	IP21	
Lifetime	5 years	
Protection Against Electric Shock	Internally powered supply	
Operating Environment	Temperature condition	5 °C - 40 °C
	Humidity condition	15 % - 90 % RH
	Atmospheric condition	70 kPa - 106 kPa
Transportation and Storage Environment	Avoid strong impact, direct impact, exposure or rain during transportation. The device shall be stored indoors at the temperature of -20 °C to 55 °C and the relative humidity of 10 % - 93 %. Atmospheric condition: 70 kPa - 106 kPa without corrosive gas and with good ventilation.	

AUTHORIZED COMPONENT

Please use the authorized adapter. Adapter is not included in the scope of delivery.

Adapter

Type: SINGOF-5U-050100
Input: 100 - 240 V - 50/60 Hz 0.16 - 0.07 A
Output: 5.0 V $\overline{\text{---}}$ 1.0 A

APPENDIX: EMC INFORMATION

Guidance and manufacturer's declaration - Electromagnetic emission

The Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Blood Pressure Monitor should assure that it is used in such an environment.

Emissions	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable	

APPENDIX: EMC INFORMATION

Guidance and manufacturer's declaration - electromagnetic Immunity									
The Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Blood Pressure Monitor should assure that it is used in such an environment.									
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)		Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385		380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27	27
	450		430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710		704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9	9
	745								
	780								
	810		800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE-Band 5	Pulse modulation 18 Hz	2	0.3	28	28
	870								
	930								

APPENDIX: EMC INFORMATION

Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)				1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
				1845	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
				1970							
				2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
				5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9	
				5500							
				5785							

APPENDIX: EMC INFORMATION

Guidance and manufacturer's declaration - Electromagnetic immunity

The Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power frequency Magnetic field IEC 61000-4-8	30A/m, 50/60Hz	30A/m, 50/60Hz
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

NOTE: UT is the a.c. mains voltage prior to application of the test level

EN

Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	30 kHz	CW	8	8
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 kHz	Pulse modulation 50 kHz	7.5	7.5

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Shenzhen AOJ Medical Technology Co., Ltd.
Share Info GmbH
Share Info Suisse GmbH
Share Info Ltd
Geratherm Medical GERMANY GmbH