Suprasorb® P

PU-Schaumverband, nicht klebend
Pansement hydrocellulaire PU, non adhésif
PU foam dressing, non-adhesive
Apósito de espuma de PU, no adhesivo
Penso de espuma de PU, não adesivo
Medicazione in schiuma di PU, non adesiva
PU-schuimverband, niet-klevend
PU-skumforbinding, ikke klæbende
PU-skumförband, icke-självhäftande
Polyuretan-skumbandasje, ikke-klebende
PU-vaahtoside, kiinnittymätön
PU pěnový obvaz, nelepící
PU penové krytie, nelepivé
Nem tapadó PU habkötszer
Opatrunek z pianki poliuretanowej, nieprzylepny
PU penasta obloga, nelepljiva
PU-pjenasti povoj, neljepiv
Αφρώδες επίθεμα PU, χωρίς κόλλα
Pansament din spumă PU, neadeziv
PU-vahthaavaside, mittekleepuv
PU putų tvarstis, nelipnus
Yapışkansız PU köpük örtü
ضمادة البولي يوريثين الرغوية، غير اللاصقة

zuschneidbar · découpable · can be cut to size · recortable · pode ser recortado · ritagliabile · kan op maat worden geknipt · kan tilskæres · kan klippas till · kan klippes · leikattava · lze přistřihnout · možné zostrihnúť na požadovanú veľkosť · méretre vágható · możliwość przycięcia · se lahko prireže · može se rezati · κόβεται · poate fi tăiat · parajaks lõigatav · galima karpyti pagal dydį · kesilebilir · يمكن قصّها

STERILE EO

0123 MD PAP

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en Instructions for Use

Product description
The Suprasorb® P PU foam dressing, non-adhesive, is a three-layer wound dressing made from a polyurethane backing layer, a polyurethane foam layer and a perforated polyurethane wound contact layer.
With regard to its fluid absorption capacity, the Suprasorb® P PU foam dressing, non-adhesive, is designed for treating moderately exuding superficial wounds. The dressing conforms to the wound bed. It maintains a moist wound healing environment, assisting the natural healing process.*

Dressing changes are gentle to the tissue as adhesion to the wound is prevented for the most part. The Suprasorb® P PU foam dressing, non-adhesive, is soft and flexible and therefore also suitable for areas of the body which are difficult to dress.

The Suprasorb® P PU foam dressing, non-adhesive, is also suitable as a secondary dressing, in combination with the following products:

- Suprasorb® A Pro,
- Suprasorb® Liquacel,
- Lomatuell® Pro.

Suprasorb® P PU foam dressing, non-adhesive, can be used as a secondary dressing for infected wounds and wounds at risk of infection as part of a regular antimicrobial treatment regimen. This product can also be used in combination with commonly used wound irrigation solutions and compression and conforming bandages.

*based on in-vitro data

Product composition
Polyurethane-film, Co-polyamide, Polyurethane-foam, Acrylic-adhesive, Polyurethane-film

Intended purpose
The PU foam dressing Suprasorb® P is intended for absorption of wound exudate and mechanical protection of the wound.

Suprasorb® P PU foam dressing, non-adhesive, is used to treat

moderately exuding,

superficial wounds.

Indications
Acute wounds:
- skin graft donor sites
- surgical wounds

Chronic wounds:
- pressure ulcers
- lower leg ulcers
- diabetic ulcers

Contraindications
Known allergy and/or hypersensitivity to any of the product components

Side effects
In individual cases skin irritations and/or allergies have occurred.

Warnings and precautions
- Suprasorb® P PU foam dressing, non-adhesive, is intended for single use and must not be resterilised.
- Do not use the Suprasorb® P PU foam dressing, non-adhesive, with oxidising substances, such as hypochlorite solutions or hydrogen peroxide.
- Do not use if package is damaged or accidentally opened.

Application instructions
For application by healthcare professionals.
When treating venous lower leg ulcers, the prescribed compression therapy, and in the case of pressure ulcers a pressure relief should be continued.

1. Prepare the wound:
Clean the wound carefully in accordance with current clinical standards. The skin surrounding the wound must be clean and dry.

2. Applying the wound dressing:
a) The size of the dressing depends on the wound area. Can be cut to size if required. The dressing should overlap the wound by 2-3 cm at least.
b) Check the sterile barrier for obvious damage prior to opening.
c) Open the packaging at the place provided for this purpose and remove the product, keeping it sterile.

wound-facing side

d) Place the dressing on the wound and secure it with a suitable secondary dressing.

3. Dressing change:
The intervals at which the Suprasorb® P PU foam dressing, non-adhesive, must be changed are determined by the treating physician based on the wound condition.
a) The Suprasorb® P PU foam dressing, non-adhesive, can remain on a carefully cleaned wound for up to 7 days.
b) The dressing should be changed when the exudate in the saturated foam layer extends beyond the size of the wound area.

Disposal
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions
In case of serious incident, contact the manufacturer and the responsible health authorities.

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