Instructions for Use

EN

Cutimed® Siltec® PLUS Heel

Soft tack silicone foam dressing with superabsorbers

Product description

Cutimed® Siltec® PLUS Heel is a silicone coated foam dressing. Cutimed® Siltec® PLUS Heel is a sterile, single-use absorbent polyurethane foam dressing for atraumatic dressing changes that contains superabsorbent stripes which absorb and lock wound exudate. The wound contact surface is a perforated silicone layer that allows adherence to the periwound skin but not to the moist wound bed or to newly formed epithelial tissue. This minimizes trauma and pain during dressing changes. The outer film is water-repellent, yet permeable to oxygen and vapor.
The foam dressing and the superabsorbent stripes are designed to absorb and to lock away excess wound fluid inside the dressing, promoting a moist wound environment and minimizing the risk of maceration. Cutimed® Siltec® PLUS Heel is highly conformable and provides protection of the wound site. Cutimed® Siltec® PLUS Heel leaves no residues in the wound. The dressing reliably retains wound exudate, even under compression. If necessary, the dressing can be cut to size. Cutimed® Siltec® PLUS Heel is anatomically shaped for the treatment of heel wounds, gently adheres to the periwound skin and stays in place while secondary fixation is applied.

Intended purpose

Cutimed® Siltec® PLUS Heel is intended for the treatment of exuding wounds with low to high exudate levels such as venous and arterial ulcers, pressure ulcers, diabetic foot ulcers, surgical incisions, skin grafts and donor sites, lacerations or abrasions.
Cutimed® Siltec® PLUS Heel may assist in the prevention of pressure ulcer as part of a comprehensive plan of continuous care, risk assessment and preventive care by a healthcare facility and healthcare professional.

Please contact your healthcare professional if you are unsure whether the product is appropriate for you.

Contraindications

There are no known contraindications.

Warnings and precautions

Cutimed® Siltec® PLUS Heel is packaged for single-use. Do not re-use or re-sterilise as there is a risk of transmission of body fluids or contaminated tissue between patients. Do not use if the pouch is already open or damaged as the sterility of the device is guaranteed only when the pouch is unopened and undamaged prior to use. Discard open or unused material. Do not use with oxidizing solutions such as hypochlorite or hydrogen peroxide. Remove the dressing prior radiation therapy. For external use only.
The wound should be inspected for signs of infection and treated according to clinical practice if required. In rare cases skin reactions (e.g. redness, itching) may occur.
If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate healthcare professional.

Directions for use

1. Clean the wound and your hands according to good clinical practice. If you are unsure how to clean the wound appropriately, please contact your healthcare professional for further instructions.
2. Make sure the surrounding skin is dry.
3. Make sure that Cutimed® Siltec® PLUS Heel covers the wound and overlaps the wound edges. (A)
4. Please ensure, that the pouch is not open or damaged. Open the pouch and take out the dressing.
5. Remove protective liner film and apply with the adherent side towards the wound. (B) (C)
6. Avoid stretching the dressing.
7. Apply additional dressing fixation.
8. Cutimed® Siltec® PLUS Heel may be used on an infected wound in addition to antimicrobial therapy under medical supervision.
9. If necessary, Cutimed® Siltec® PLUS Heel can be used in combination with other dressings e.g. Cutimed® Sorbact®/Sorbact® (in case of infected wounds) or hydrogel e.g. Cutimed® Gel (for dry, necrotic wounds).
10. Cutimed® Siltec® PLUS Heel can remain in place for up to 7 days, depending on wound conditions and saturation of wound dressing.
 It is recommended to change the dressing every 24 hours initially, moving to less frequent changes as necessary.
11. Carefully remove any additional dressing fixation and gently remove Cutimed® Siltec® PLUS Heel from the wound and discard.
12. If required, cleanse the wound prior to application of a new dressing.

Storage and Disposal

Store dressing away from direct sunlight at ambient temperature and humidity. Keep out of the reach of children.
Discard used product according to local regulations, e.g. via regular household or clinical waste.

Notice

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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