Instructions for Use

Cutimed® Gelling Fiber

(en)

Product Description
Moisture maintaining gelling fiber dressing.
Cutimed® Gelling Fiber is a soft, hydrophilic, sterile, nonwoven, gel forming pad or ribbon dressing composed of sodium carboxymethylcellulose and enhanced fibres for strengthening. Cutimed® Gelling Fiber is a conformable and highly absorbent dressing. Cutimed® Gelling Fiber absorbs wound fluid and transforms into a soft gel, which maintains a moist environment to support the body's healing process. Cutimed® Gelling Fiber helps in the removal of nonviable tissue from the wound (autolytic debridement) without damaging newly formed tissue.

Intended purpose
Cutimed® Gelling Fiber is intended for the treatment of moderately to highly exuding wounds healing by secondary intention, such as: leg ulcers, pressure ulcers (stage II-IV) and diabetic ulcers, surgical wounds (e.g., post-operative, wounds left to heal by secondary intent and donor sites), partial thickness burns, exudate absorption in oncology wounds, traumatic wounds (e.g., abrasions and lacerations). The dressing can be used for dry wounds, but care should be taken when removing the dressing. Cutimed® Gelling Fiber can be used under compression.
Please contact your healthcare professional if you are unsure whether the product is appropriate for you.

Contraindications
There are no known contraindications.

Warnings and Precautions
Do not use if the pouch is already open or damaged as the sterility of the dressing is guaranteed only when the pouch is unopened and undamaged prior to use. Cutimed® Gelling Fiber is for single use only. Re-use may lead to increased risk of infection or cross contamination. Do not re-sterilize. Cutimed® Gelling Fiber should be used under medical supervision. Please always follow the instructions of any product that is used in combination. If the treated condition deteriorates (e.g., development of infection), fails to improve or if a side effect is observed, consult a physician or an appropriate health care professional. Because Cutimed® Gelling Fiber creates an environment that favors new blood vessels' growth, delicate newly formed blood vessels may occasionally produce a blood-stained wound fluid. Cutimed® Gelling Fiber is not intended for use as a surgical sponge. Cutimed® Gelling Fiber is not intended for use within internal body cavities or within closed wounds. The wound should be inspected for signs of infection and treated according to clinical practice. If required, please contact your healthcare professional.

Directions for Use
Application and removal of the dressing should be in accordance with good clinical practice.

Preparation
1. Before applying the dressing disinfect your hands and use gloves. Prepare treatment environment and items for an aseptic dressing application (sterile tweezers and scissors). If necessary, debride (e.g. with Cutimed® DebriClean) according to the instructions for use, irrigate the wound and remove leftover cleansing solution afterwards. Select a suitable size of the dressing.
2. Remove Cutimed® Gelling Fiber from pouch wearing sterile gloves or using sterile tweezers (**Fig. A**)
3. Cut to size if necessary.
4. Both sides of the product can be applied to the wound.

Application
Exuding and dry wounds (excluding partial thickness burns)
- Cutimed® Gelling Fiber should overlap the wound by at least 1 cm, as the dressing will shrink when it absorbs exudate. (**Fig. B1**)
- Apply the dressing to the exuding wound and cover it with a secondary moisture balancing dressing such as Cutimed® Siltec® (B / Sacrum) or Cutimed® Sorbion® Border. For low exuding wounds a secondary fixation such as Hypafix® transparent or Fixomull® transparent can be used. (**Fig. C**)
- For dry wounds apply the dressing onto the wound. Wet the dressing with sterile water or saline solution over the wound area and cover it with a moisture balancing secondary fixation such as Hypafix® transparent or Fixomull® transparent to avoid drying out of the dressing and subsequent dressing adherence to the wound. (**Fig. C**)

Deep wounds and cavities
- Dry the surrounding skin thoroughly.
- Loosely pack ribbon or pad into the wound to allow swelling. When using Cutimed® Gelling Fiber ribbon, cut to appropriate length leaving a small overhang of 2 - 3 cm outside the wound for easy removal. (**Fig. B2**)
- Cover with a secondary moisture balancing dressing such as Cutimed® Siltec® (B / Sacrum) or Cutimed® Sorbion® Border. (**Fig. C**)

Partial thickness burns and donor sites
- For larger wound areas, select one or more dressings that the dressing overlaps the wound site by at least 5 cm as the dressing will shrink when absorbing exudate. (**Fig. B1**)
- Cutimed® Gelling Fiber should be covered with a sterile, loosely woven gauze pad and appropriately secured. (**Fig. C**)
- Change the cover regularly and inspect Cutimed® Gelling Fiber while it remains in place on the wound.
- Adherence to the wound bed of Cutimed® Gelling Fiber is a desired characteristic. Nonadherence of areas of the dressing may indicate deepening of the wound or infection. Areas of the dressing may be cut away to enable inspection of the wound. The exposed areas should then be treated appropriately.
- As the wound re-epithelializes, Cutimed® Gelling Fiber will detach or can easily be removed. Cutimed® Gelling Fiber may be used in combination with a Sorbact-technology product as a wound contact layer, if indicated.

Dressing change
All wounds should be inspected regularly. Remove Cutimed® Gelling Fiber when clinically indicated (e.g. excessive bleeding, leakage or suspicion of infection). Cutimed® Gelling Fiber can be left in place for up to 7 days, except in burns and donor sites where it can be left in place for up to 14 days. Partial thickness burns and donor sites: In the post wound period (up to 4 days), high volumes of wound exudate may require frequent changes of the dressing.

Removal
If removing Cutimed® Gelling Fiber is difficult, it should be fully saturated with sterile saline or water and removed carefully and slowly. Cleanse the wound after removal of the dressing to avoid residues with an appropriate sterile wound cleansing solution and remove the leftovers.

Supportive measures
Appropriate supportive measures should be taken where indicated e.g. for venous leg ulcers use of graduated compression bandaging (e.g. JOBST® Ulcer Care)

Size table

cm	inch	cm	inch
5 x 5 cm	2 x 2 in.	4 x 20 cm	3/2 x 8 in.
10 x 10 cm	4 x 4 in.	4 x 30 cm	3/2 x 12 in.
15 x 15 cm	6 x 6 in.	2 x 45 cm	3/4 x 18 in.
4 x 10 cm	3/2 x 4 in.	1 x 45 cm	2/5 x 18 in.

Storage and Disposal
Store dressing away from direct sunlight at ambient temperature and humidity. Discard any unused portion of Cutimed® Gelling Fiber or used product according to local regulations, e.g. via regular household or clinical waste.

Notice
Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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